

Daily MAINTENANCE Check Log
 (for the **Defibtech DDU-100** Semi Automatic External Defibrillator)

AED Location: _____ AED Director: _____
 AED Serial No.: _____

Month:		Is the Active Status Indicator flashing green?		
Day	Initial	Yes	No	If No, see manual for repair options and enter the date (if any) the AED is returned for service
1		<input type="radio"/>	<input type="radio"/>	
2		<input type="radio"/>	<input type="radio"/>	
3		<input type="radio"/>	<input type="radio"/>	
4		<input type="radio"/>	<input type="radio"/>	
5		<input type="radio"/>	<input type="radio"/>	
6		<input type="radio"/>	<input type="radio"/>	
7		<input type="radio"/>	<input type="radio"/>	
8		<input type="radio"/>	<input type="radio"/>	
9		<input type="radio"/>	<input type="radio"/>	
10		<input type="radio"/>	<input type="radio"/>	
11		<input type="radio"/>	<input type="radio"/>	
12		<input type="radio"/>	<input type="radio"/>	
13		<input type="radio"/>	<input type="radio"/>	
14		<input type="radio"/>	<input type="radio"/>	
15		<input type="radio"/>	<input type="radio"/>	
16		<input type="radio"/>	<input type="radio"/>	
17		<input type="radio"/>	<input type="radio"/>	
18		<input type="radio"/>	<input type="radio"/>	
19		<input type="radio"/>	<input type="radio"/>	
20		<input type="radio"/>	<input type="radio"/>	
21		<input type="radio"/>	<input type="radio"/>	
22		<input type="radio"/>	<input type="radio"/>	
23		<input type="radio"/>	<input type="radio"/>	
24		<input type="radio"/>	<input type="radio"/>	
25		<input type="radio"/>	<input type="radio"/>	
26		<input type="radio"/>	<input type="radio"/>	
27		<input type="radio"/>	<input type="radio"/>	
28		<input type="radio"/>	<input type="radio"/>	
29		<input type="radio"/>	<input type="radio"/>	
30		<input type="radio"/>	<input type="radio"/>	
31		<input type="radio"/>	<input type="radio"/>	

Check the Status Indicator to ensure that it is flashing **GREEN**. When the indicator is **GREEN**, the AED is ready for a rescue.